

- If you are viewing this course as a recorded course after the live webinar, you can use the scroll bar at the bottom of the player window to pause and navigate the course.
- This handout is for reference only. Non-essential images have been removed for your convenience. Any links included in the handout are current at the time of the live webinar, but are subject to change and may not be current at a later date.

No part of the materials available through the continued.com site may be copied, photocopied, reproduced, translated or reduced to any electronic medium or machine-readable form, in whole or in part, without prior written consent of continued.com, LLC. Any other reproduction in any form without such written permission is prohibited. All materials contained on this site are protected by United States copyright law and may not be reproduced, distributed, transmitted, displayed, published or broadcast without the prior written permission of continued.com, LLC. Users must not access or use for any commercial purposes any part of the site or any services or materials available through the site.

Technical issues with the Recording?

- Clear browser cache using [these instructions](#)
- Switch to another browser
- Use a hardwired Internet connection
- Restart your computer/device

Still having issues?

- Call 866-782-9924 (M-F, 8 AM-8 PM ET)
- Email customerservice@OccupationalTherapy.com

Wheelchair Mobility: Power Add On Options And Use For Manual Wheelchairs

Michelle L. Lange, OTR/L, ABDA, ATP/SMS



Quickie Xtender

Learning Outcomes

The participant will be able to:

1. Define Alternate Drive Mechanisms and provide examples of 3 types.
2. Describe clinical indicators for Alternative Drive Mechanisms.
3. Describe clinical indicators for PAPAWs.

What we will be covering:

- Alternate Drive Mechanisms
- Providing assist for self-propulsion of a manual wheelchair
 - Lever drives
 - Arm-crank drives
 - Geared hubs
 - Reverse pushrim drives
 - Push rim-activated power assist drives

3

SMS Series

- This is part of a series of webinars designed to prepare the participant for the Seating and Mobility Specialist examination
- And... develop more advanced seating and wheeled mobility skills



4

Mobility Hierarchy

- Dependent Manual Wheelchairs
- ↓
- Independent Manual Wheelchairs
 - Designed for self-propulsion
- ↓
- **Alternate Drive Mechanisms**
- ↓
- Power Wheelchairs

5

Alternate Drive Mechanisms

- Goal:
 - Improved mechanical efficiency
 - Maintain transportability, versatility, and aesthetics of a manual wheelchair
 - DiGiovine & Berner, 2018

NuDrive Air



6

CONTINUED

Who is this Appropriate For?

- It is well documented that people who self-propel a manual wheelchair, even an appropriately configured ultralightweight MWC, are prone to:
 - Repetitive stress injuries
 - Shoulders, wrists
 - Pain
- This can worsen to the point that mobility is compromised.

7

CONTINUED

Who is this Appropriate For?

- A person who can self-propel, but:
 1. Has upper extremity injury and/or pain
 2. Has decreased cardiopulmonary function
 3. Has decreased strength and endurance

Quickie Q7



8

CONTINUED

Who is this Appropriate For?

- A person who can self-propel, but:
 4. Is less efficient in mobility secondary to aging with a disability or aging in general
 5. Is a long-term MWC user who may have lost efficiency
 6. Is a MWC user who cannot manage slopes, varied terrain and/or long-distances without assistance

Quickie Q7



9

Non-powered options

- Geared systems
- Lever systems
- Arm-crank systems

10

Geared Systems

11

Geared Systems

- Gear system between the handrim and the wheel
 - Geared hubs
 - Reverse pushrim drives




Intelliwheels
Easy push
wheel

12

Geared Systems


- How does it work?
- Standard MWC handrim has a 1:1 ratio
 - 1 push = set distance of movement



13

Geared Systems

- Geared systems
 - Typically 1:1 ratio for standard propulsion
 - 1:2 ratio for propulsion on **inclines** or **high rolling resistance surfaces**
 - i.e. carpet
 - Reduces UE effort on all ramps
 - Reduces abdominal activity on steep ramps



14

Geared Systems

- Who is this appropriate for?
 - Users who could not traverse certain environments without this system
 - Users with reduced UE strength

15

Geared Systems

- Reverse propulsion
- Pulling back on the wheel moves the MWC forward



RoWheels
Revolution

16

Geared Systems

- Pulling action uses muscles in your upper back and shoulders
 - Larger muscles
 - Encourages better posture
 - Reduced risk of shoulder-related injuries
 - Minimizes impact on handrim



RoWheels
Revolution

17

Geared Systems

- The user must recognize the need to change gears
- Adds weight
 - Most weight is by the hub, less impact on propulsion
- Minimizes roll-back on slopes
- Reduces force required to brake when going downhill

18

Lever Systems

19

Lever Systems



- Lever arm increases torque compared with traditional handrim
- Typically levers are placed between frame and wheel
- Does not add very much to width of chair
- Gears may also be built-in
- Handle is easier to grasp than the handrim

20

CONTINUED

Lever Systems

- Steering
 - Differential motion
 - Pull back on one side, push forward on the other for a turn
 - Dedicated steering mechanism
 - Turn the handle of the lever



GRIT
Freedom
Chair

21

CONTINUED

Lever Systems

- More mechanically efficient than standard handrims
- Speed higher
- Heart rate and O_2 consumption less
- Less demanding, less strenuous, more efficient



NuDrive Air

22

CONTINUED

Armcrank Systems

23

Armcrank Systems

- Often used in handcycles
- User may be upright or recumbent
- Gears similar to a bicycle
- Cons:
 - Bulky
 - Reduced maneuverability
 - Aesthetics



24

Armcrank Systems

- Most systems are synchronous
 - Arms turn together – more efficient
- Continuous movement is less straining and more efficient than handrims
- Longer distances and higher speeds
- Recumbent position reduces air resistance

Invacare Top
End XLT



25

Armcrank Systems

- Typically used for recreational purposes
- Although efficient, these are not practical for everyday use in a variety of environments

26

Justification

- How are these systems justified?
 - Allows a MWC user to be maintain independence in mobility without having to switch to a PWC
 - Increases efficiency while decreasing risk of injury and pain
 - Cost of these systems is typically less than a PWC



Quantum Edge 3

27

Powered Options

28

CONTINUED

Powered options

- Pushrim-Activated Power Assist (PAPAW)

- Motor in hub of the wheel
- Motor external to the frame

Quickie
Xtender



Max Mobility
SmartDrive



29

CONTINUED

PAPAWs

- How does it work?

- Hybrid system
- Measures force imparted on pushrim and amplifies force through electromechanical system
- Amplifies the propulsion phase (faster and increased force)
- Extends the recovery phase (time before another stroke required)

Quickie
Xtender
battery



30

CONTINUED

CONTINUED

PAPAWs

- So what exactly does that mean?
 - A typical push on the handrim results in more force to propel the MWC further per stroke and with more force up slopes and over varied terrain
 - Cool!

31

CONTINUED

PAPAWs

- Modes of operation
 - Intermittent
 - Continuous mode

32

CONTINUED

CONTINUED

PAPAWs

- Intermittent mode
 - Motor engages only in response to force on handrim
 - Only used when needed, may be safer for some users

Frank Mobility
E-Motion



33

CONTINUED

PAPAWs

- Continuous mode
 - Motor is turned on
 - Via Switch or BT device (i.e. wristband)
 - Handrims only used for course corrections
 - Great for long distances, requires minimal effort
 - Does require attention

Max Mobility
SmartDrive



34

CONTINUED

PAPAWs

- Translating handrim forces into direction of movement
 - Some systems react differently to forces on each handrim. User may push with uneven forces. Most similar to typical MWC propulsion.
 - Some systems read forces on both handrims together for smoother propulsion. Smarter electronics.



Frank
Mobility
E-motion

35

PAPAWs

- Uses
 - Client who is not an efficient MWC user, but is reluctant to move to a PWC
 - Injury prevention
 - Energy conservation
 - Increases distance for some users
 - Increases access/control on slopes and varied terrain

36

PAPAWs

- Drawbacks
 - Weight
 - Cost
 - Reduces ability to do wheelies
 - May impact accessibility



Max Mobility
SmartDrive

37

So if these increase efficiency that much...

- Why doesn't everyone use them?
- Good question!
- Research has shown that all these strategies basically compensate for MWC design issues that lead to injury and pain that worsens over time.

38

So if these increase efficiency that much...

- Most users are not aware of these technologies and are often resistant to more technology than they need.
- Most funding sources are not aware of these technologies and are reluctant to pay for more than they believe is necessary.
- Future – hopefully improved design and awareness through education

Case Study

Case Study

- DeShawn
- Transverse Myelitis
- 13 years old
- Middle School



Quickie 2

41

DeShawn


- DeShawn was using an ultralightweight MWC with customized seating.
- He was able to propel throughout his school, though this was slow and took a lot of energy, especially for longer distances.
 - Moving from class to class
 - Led to fatigue that was impacting school performance

42

CONTINUED

DeShawn

- He also could not slow his speed when going down a ramp.
- He did not have adequate strength to propel up a ramp.
- As a result, he was dependent for these activities.



43

CONTINUED

DeShawn

- Power Wheelchair?
- We considered a power wheelchair, but DeShawn was very reluctant to move to a power wheelchair.
- The family did not have an accessible vehicle, and they needed something they could fold.

Quickie QM-7



44

CONTINUED

DeShawn

- We arranged for a trial with some PAPAWs
- We could remove his rear wheels so that he could try the PAPAWs on his own MWC.
- He liked them!



Quickie
Xtender

45

DeShawn

- Funding
 - We were concerned about funding, as these were costly
- We measured how long it took DeShawn to move between specific locations before and after the PAPAWs.
- We documented situations he was unable to manage in the current MWC.
 - i.e. up a ramp
- We then documented if he could perform that task with the PAPAWs.

46

DeShawn

- Results
 - DeShawn was independent in all of his mobility, including varied terrains and slopes.
 - He no longer had to depend on others or risk becoming 'stuck'.
 - He was not as fatigued and so was better able to participate at school.

49

Take Home Message:

- Many people using a MWC are inefficient in self-propulsion.
- Many people who self-propel develop injuries and pain.
- Options exist which close the gap between MWCs and PWCs.

50

References:

- DiGiovine, C. & Berner, T. (2018). In Seating and Wheeled Mobility: a clinical resource guide, eds M. Lange & J. Minkel, Slack, Thorofare, NJ.

51

Questions?

- Email

52

CONTINUED

Thank You!

CONTINUED

Contact Information

- Michelle L. Lange, OTR/L, ABDA, ATP/SMS
- MichelleLange1@outlook.com
- www.atilange.com

CONTINUED