

- If you are viewing this course as a recorded course after the live webinar, you can use the scroll bar at the bottom of the player window to pause and navigate the course.
- This handout is for reference only. It may not include content identical to the PowerPoint. Any links included in the handout are current at the time of the live webinar, but are subject to change and may not be current at a later date.

© continued.com, LLC 2017. No part of the materials available through the continued.com site may be copied, photocopied, reproduced, translated or reduced to any electronic medium or machine-readable form, in whole or in part, without prior written consent of continued.com, LLC. Any other reproduction in any form without the permission of continued.com, LLC is prohibited. All materials contained on this site are protected by United States copyright law and may not be reproduced, distributed, transmitted, displayed, published or broadcast without the prior written permission of continued.com, LLC. Users must not access or use for any commercial purposes any part of the site or any services or materials available through the site.

Dynamic Seating

Michelle L. Lange, OTR/L, ABDA, ATP/SMS



Learning Objectives

1. The participant will be able to list 3 clinical indicators for a dynamic movement of the pelvis.
2. The participant will be able to list 3 clinical indicators for dynamic movement of the lower extremities.
3. The participant will be able to list 3 clinical indicators for dynamic movement of the neck.

What we are covering today:

- What is Dynamic Seating?
- Pelvis
- Lower Extremities
- Head



Dynamic Seating – a definition

- Dynamic Seating is movement which occurs within the seat and/or wheelchair frame in response to force from the client. Dynamic components absorb force which in turn assists the client back to a starting position.



Dynamic Seating – Goals

- Primary Goals:

1. To allow movement
2. To diffuse force
3. To protect the client, seating system, mounting hardware, and mobility base
4. To improve postural control



Goal #1: To Allow Movement

- What are the benefits of movement in the wheelchair?

- To increase sitting tolerance and compliance
- To provide vestibular input
- To increase alertness
- To decrease agitation
- To increase function
- To provide active range of motion
- *video



Goal #2: To Diffuse Force

- By diffusing force, we achieve these goals:

- *To reduce active extension
- To reduce energy exertion

- Which in turn, may also help:

- To increase sitting tolerance and compliance
- To decrease agitation
- To increase function



Goal #3: To Protect

- To protect the client

- If the client is exerting enough force to break components, injury is very possible
- Micro-concussions
- Other injuries

- To protect the seating system, mounting hardware, mobility base frame



Goal #4: Postural Control

- By providing movement against light force, strength may build
- This can lead to improved trunk and head control
- **Warning:** the primary purpose of dynamic seating is not therapeutic



Dynamic Seating: the pelvis

- Allowing movement at the pelvis has advantages and disadvantages



Pelvic Dynamic Seating: advantages

- If pelvic movement is blocked, this force can be transferred to other body areas, resulting in increased extension. Providing movement at the pelvis reduces overall extension.
- Movement of the pelvis shifts weight which provides pressure relief and comfort



Pelvic Dynamic Seating: disadvantages

- Movement may open seat to back angle which could result in a posterior pelvic tilt
 - This may be acceptable is the pelvis returns to neutral upon return to upright



Pelvic Dynamic Seating: disadvantages

- Allowing movement of the pelvis can lead to **assumption of a destructive posture**
- Allowing movement of the pelvis into posterior pelvic tilt can lead to **increased extension and spasms**
- The client may not be able to **return to a neutral position**



Pelvic Dynamic Stability: product options

- Integrated systems
- Dynamic Backs
 - Miller's Adaptive Equipment
 - Otto Bock
 - Seating Dynamics
 - Stealth
 - Sunrise Medical

Integrated Systems

- These systems are a complete seat and wheelchair, generally, and not retrofittable to other bases
- Multiple dynamic movements work together
- Most are only available outside of the United States

Integrated Systems available outside the US

- Adacta Klim
- Aktivline
- Netti



continued

Elevation

- Elevation Ultra Lightweight wheelchair
- PDG
- Seat can rise at rear, 10 degrees in height
- 30 degrees back rest angle adjustment



continued

Leggero Dyno

- Activator Dynamic Seating Component
 - Active knee range 15 degrees
 - Active seat to back angle 35 degrees
- Folds
- *video



continued

Kinetic Innovative Seating System

- KiSS for Wheelchairs, in USA (not yet available)
- Allows constant articulated motion
- Can fit on most manual wheelchair frames
- Seat Back and Seat Base can be installed separately or together



continued

Dynamic Backs

- Movement occurs only at the back
- Can often be combined with other dynamic options to provide movement in other areas

continued

Miller's Dynamic Backrest Interface

- Extends at level of biangular back
- *video

continued

Miller's Dynamic Back

- *video

CONTINUED

Otto Bock Dynamic Back

- Dynamic Extension Relief Back
- 4 sets of springs between back and mounting hardware
- Captures movement back and rotation
- Springs are shrouded



CONTINUED

Seating Dynamics Dynamic Back

- Seating Dynamics
 - Dynamic rocker back
 - Resistance is adjustable through a set of elastomers



CONTINUED

Seating Dynamics Dynamic Back

- Seating Dynamics
 - Dynamic rocker back
 - *video
 - *video

CONTINUED

Stealth Products Dynamic Back

- Stealth Dynamic Backrest Mounting Hardware
 - Encased to protect mechanism



CONTINUED

Stealth Products Dynamic Back

- Tarta back
- Off the shelf, can customize
- Similar to the Ortho Flex back
 - Made in Italy, but available in the USA
- Goal: to assist movement to improve function
- *video



CONTINUED

Sunrise Medical Dynamic Back

- Mono Back or Dual Cane
- Available on Quickie manual wheelchairs
- Locks out
- Dynamic option



continued

Questions?

continued

Dynamic Seating: the lower extremities



continued

LE Dynamic Seating: advantages

- Many clients will not tolerate having their feet restrained
- Stability is often required at the feet, however, to improve function
- Dynamic seating may improve tolerance and compliance, while providing function
- Limiting lower extremity movement may protect the feet from injury



continued

LE Dynamic Seating: disadvantages

- Restricting the feet in any way will prevent independent transfers
- Some clients will continue to fight any restraint of the feet



continued

Lower Extremity Dynamic Seating:
product options

- Miller's Adaptive Equipment
- Seating Dynamics

continued

Miller's Dynamic Footrests

- Dynamic Footrest Gas Spring
 - Extends downward 2"
- Dynamic Footrest Coil
 - Allows some rotation at the footplate



CONTINUED

Miller's Dynamic Footrests

- Dynamic, Articulating Footrest Hanger
- *video

CONTINUED

Seating Dynamics Dynamic Footrests

- Dynamic footrest
- Telescopes 1 1/2"
 - Use alone for clients who may otherwise lose the position of their pelvis
 - Tight hamstrings
- Optional knee extension 30 degrees
- Optional dynamic dorsi/plantar flexion, 17 degrees each direction
- Resistance is adjustable through springs



CONTINUED

Seating Dynamics Dynamic Footrests

- Dynamic footrest
- *video
- *video

CONTINUED

Questions?

CONTINUED

Dynamic Seating: the head



CONTINUED

Head Dynamic Seating: advantages

- Providing some movement can:
 - Absorb force and protect the neck and brain
 - Reduce breakage of head support mounting hardware
 - Reduce loss of alignment of head support
 - Diffuse force



Head Dynamic Seating: disadvantages

- Movement can lead to postural insecurity
- Excessive movement can trigger reflexive response
 - Moro
 - Tonic neck



Head Dynamic Seating: product options

- Metalcraft
- Miller's Adaptive Equipment
- Otto Bock
- Seating Dynamics
- Stealth Products
- Symmetric Designs

Dynamic Headrest Options

- Metalcraft
 - Bi-directional
 - Prototype
 - *video

Dynamic Headrest Options

- Miller's Dynamic Headrest Interface
- *video

Dynamic Headrest Options

- Miller's Dynamic Headrest Horizontal Adjustment Bar
- *video

Dynamic Headrest Options

- Otto Bock Dynamic Rock-n-Lock Headrest Bra
- Spring loaded mechanism, 1 1/2" travel
- Shrouded to protect hands and hair



Dynamic Headrest Options

- Seating Dynamics Dynamic Headrest
 - Single Axis moves along midline or the Y Axis, 8 degrees
 - Resistance can be changed using different elastomers
 - Multi-Axis moves in both X and Y Axis and anywhere in between
 - Capturing posterior and rotational movements
- *video



Dynamic Posterior Head Supports

- Stealth Tone Deflector
 - 10 degrees any direction
 - Works well for clients who do not tolerate a larger degree of movement
 - Protects hardware
 - Absorb and Avert!
- *video



TD-100

Dynamic Headrest Options

- Symmetric Designs Axion Rotary Interface
 - Friction knob



Dynamic Headrest Options

- Whitmyer Flex Interface Bracket
 - Provides movement upward and then back to neutral



Movement and Dynamic Components at the head

- Phillip doesn't have dynamic head components...yet
 - Note bald spot!
- *video

Questions?

Combination Approach

- Remember, these components can be used in combination
- *video

Daniel

- Daniel
- Cerebral palsy
- Age 9
- Manual tilt in space wheelchair
- Linear seating system



Daniel

- The Problem:
- Daniel is extremely strong. He routinely breaks seating components, has dislocated both elbows and has injured his knees from strong extension
- He has a Baclofen pump, but cannot tolerate increased doses due to seizures

Daniel

- Kid Rock
- Daniel trialed a Kid Rock for 2 weeks. He liked this system and could easily engage the springs
- The spring tension in the back was inadequate to consistently return to upright
 - Stronger springs



Daniel

- Daniel did very well in a Kid Rock 2
- He eventually moved into an Aspen Seating Orthosis (2 piece) as he was beginning to develop a scoliosis



Daniel

- As he grew, he was out of proper alignment with the pivot points and the system was no longer meeting his needs
- The Kid Rock 3 fits him, but is so large that the family returned it
- He has a new tilt in space MWC with Seating Dynamic components
- This has met his needs very well!



CONTINUED

Questions?

CONTINUED

Take home message:

- Dynamic Seating can either allow movement of the client within the seating system or provide movement of the seating component and/or frame
- Dynamic Seating can protect the seat and frame from damage by diffusing force
- Dynamic Seating can protect the client from undue forces and reduce tone and posturing by diffusing force
- Dynamic Seating can provide active movement

CONTINUED

Thanks!

CONTINUED

Contact Information:

Michelle Lange
MichelleLange1@outlook.com
www.atilange.com