

- If you are viewing this course as a recorded course after the live webinar, you can use the scroll bar at the bottom of the player window to pause and navigate the course.
- This handout is for reference only. Non-essential images have been removed for your convenience. Any links included in the handout are current at the time of the live webinar, but are subject to change and may not be current at a later date.

No part of the materials available through the continued.com site may be copied, photocopied, reproduced, translated or reduced to any electronic medium or machine-readable form, in whole or in part, without prior written consent of continued.com, LLC. Any other reproduction in any form without such written permission is prohibited. All materials contained on this site are protected by United States copyright law and may not be reproduced, distributed, transmitted, displayed, published or broadcast without the prior written permission of continued.com, LLC. Users must not access or use for any commercial purposes any part of the site or any services or materials available through the site.

Technical issues with the Recording?

- Clear browser cache using [these instructions](#)
- Switch to another browser
- Use a hardwired Internet connection
- Restart your computer/device

Still having issues?

- Call 866-782-9924 (M-F, 8 AM-8 PM ET)
- Email customerservice@OccupationalTherapy.com

continued

Successful Static Splinting

Hand Based Splint Fabrication, Part 3



Dr. Kirsten Davin
OTD, OTR/L, ATP, SMS

1

continued

Copyright/Permissions

- All images, charts, graphs, photos and other assets in this PowerPoint are owned by Dr. Kirsten Davin and/or she has secured written permission to publish them on continued.com.
- A special thank you to Performance Health, who provided images, charts and graphs for our educational purposes.

2

continued

Disclosures

- **Presenter Disclosure:**
 - Financial: Kirsten Davin has received an honorarium for presenting this course.
 - Non-financial: Kirsten Davin has no relevant non-financial relationships to disclose.
- **Content Disclosure:** This learning event does not focus exclusively on any specific product or service.
- **Sponsor Disclosure:** This course is presented by OccupationalTherapy.com.

Learning Outcomes

After this course, participants will be able to:

- Identify three (3) characteristics of thermoplastic materials which are applicable to hand-based splint fabrication.
- List the steps involved in the fabrication of a MP joint blocking splint and a finger gutter splint, including wear schedule.
- Describe a client's clinical presentation to determine the appropriateness of a hand-based splint application, and determine the most effective type of orthotic to implement using a case study.

continued

Materials Selection

Characteristics Beneficial to Hand Splint Fabrication

- Drapability
 - Ability of a thermoplastic to drape over a surface via the use of gravity
- Conformity
 - Ability of a thermoplastic to conform to the surface
- Memory
 - Ability of a thermoplastic to return to its original shape when reheated
- Perforation
 - Holes in the thermoplastic offer breathability

Q1

5

continued

Materials Selection

Thermoplastics Beneficial to Hand Splint Fabrication

- Aquaplast – T
 - Often preferred by pediatric clients due to color offerings and ability to accommodate
 - Multiple thickness and perforation options
 - Thermoplastic with the most ‘memory’



6

continued

Materials Selection

Thermoplastics Beneficial to Hand Splint Fabrication

- Polyform
 - The most conforming thermoplastic currently available
 - Often used for very detailed and/or small anatomical surfaces



Q2/Q3 7

Overview of Splint Fabrication

- Determine - splint and pattern
- Decide - material and how to select
- Prepare - material
- Apply - patient
- Finish - trimming, edge finishing, strapping
- Assess - fit, make necessary modifications
- Instruct - care and wearing schedule

Q4 8

continued

Splint Fabrication Goals

Serial Static

- Maintain arches
- Contour to skin
- Maintain motion (if indicated)
- Permit balanced function of unaffected muscles
- Allow maximal mobility with optimal stability
- Allows for freedom of digits
- Minimal stretch for a longer period rather than quick correction
- Larger surface area to distribute pressure following the normal contours of the hand and arm

9

continued

Hand Based Thumb Spica Splint

Design Options



10


continued

continued

MP Blockina Splint

Design Options

- Often used to treat trigger finger
- Blocks the MP joint from performing flexion, while allowing free flexion of the IP joint of the digit.



Q5/Q6


11

continued

Finger Gutter Splint

Design Options

- Often used for protection of an injured finger
- Blocks the PIP and DIP joints to keep the finger immobile and prevent flexion
- Variations in splint design



Q7

12

continued

Prefabricated Finger Splints

Oval-8 Splint

- May also be used for trigger finger to rest the tendon
- Also frequently prescribed for mallet finger
 - Occurs when the distal IP flexes and requires assistance with extension



Q8

13

Prefabricated Finger Splints

Spring Finger Extension Splint

- Designed to treat a variety of diagnoses
 - Indicated for the involvement of PIP tightness
 - Indicated for PIP joint limitations of 45 degrees or less
- Dorsal pad distributes pressure evenly across the PIP





14

Prefabricated Finger Splints

Sizing and Placing a Spring Finger Extension Splint

- Measure from the distal palmar crease to the DIP crease of the palm.
- When placing the finger extension splint, do not extend the PIP in question immediately
- Measure to determine

Spring Finger Model #501			Spring PIP Finger Model #504		
SIZE	MEASUREMENT	PRODUCT #	SIZE	MEASUREMENT	PRODUCT #
X-Small	2.25"	501AA	X-Small	1.75"	504AA
Small	2.625"	501A	Small	1.875"	504A
Medium	3"	501B	Medium	2"	504B
Large	3.5"	501C	Large	2.25"	504C
X-Large	4"	501D	X-Large	2.625"	504D




Q9 15

Prefabricated Finger Splints

Static Progressive Finger Extension Splint

- Designed to offer low, prolonged stretch to the PIP joint to encourage extension
 - Often able to be adjusted by a medically compliant client under guidance from the therapist/MD.



Q10 16

continued

Questions? Thoughts or Feedback? On-Site Course Requests?

Contact Dr. Kirsten Davin anytime!

- Phone: 217.414.2585
- E-mail: Kirs10k@aol.com
- LinkedIn: <https://www.linkedin.com/in/dr-kirsten-davin-a459a274>

17

continued

References

- AOTA (n.d.). The unique role of occupational therapy in rehabilitation of the hand. Retrieved from <https://www.aota.org/About-Occupational-Therapy/Professionals/RDP/hand-therapy.aspx>
- AOTA Fact Sheet (n.d.). Occupational therapy's role with rehabilitation of the hand. Retrieved from <https://www.aota.org/~media/Corporate/Files/AboutOT/Professionals/WhatsOT/RDP/Facts/Hand%20Therapy%20fact%20sheet.pdf>
- Baksa G., Mandl P., Benis S., Patonay L., Balint G.P., Balint P.V. (2018) Gross Anatomy of the Human Hand. In: Balint P., Mandl P. (eds) Ultrasonography of the Hand in Rheumatology. Springer, Cham. 15-41.
- Cha, Y.J., (2018). Changes in the pressure distribution by wrist angle and hand position in a wrist splint. Hand Surgery & Rehabilitation. 37(1): 38-42.

18

continued

References

- Hadidi M, Mohammad-Reza H, Kheradmand AA, Farzan M, Monsef H, et al. (2016) Is Cock-up Splint the Right Choice for All of the Carpal Tunnel Syndrome Patients? A Case Report. *Clin Med Rev Case Rep* 3:083
- Imran, K., et al. (2017). A Brief Overview of Shape Memory Effect in Thermoplastic Polymers. *Smart Polymer Nanocomposites. Springer Series on Polymer and Composite Materials*. https://doi.org/10.1007/978-3-319-50424-7_10
- Mohammad K., et al. Thermoplastic Sheet for Orthoses, A Review of Literature (2017). *EC Orthopaedics*. 189-193.
- Richards, T., Clement, R., Russell, I., & Newington, D. (2018). Acute hand injury splinting - the good, the bad and the ugly. *Annals of the Royal College of Surgeons of England*, 100(2), 92–96. <https://doi.org/10.1308/rcsann.2017.0195>

19

References

- Roll, S. C., & Hardison, M. E. (2017). Effectiveness of occupational therapy interventions for adults with musculoskeletal conditions of the forearm, wrist, and hand: A systematic review. *American Journal of Occupational Therapy*, 71, 7101180010. <https://doi.org/10.5014/ajot.2017.023234>
- Wormald, J., Millar, C., & Sivakumar, B. (2018). An Eco-conscious Alternative Material for Splinting in Hand Surgery. *Techniques in Hand & Upper Extremity Surgery*. 22(4); 156-157.

20